

Adopts UV laser printing technology, which can clearly print the information that users need on the surface of the cassette box, with a built-in touch screen, convenient operation interface, and powerful functions. KF-UE-006 is small and stable, with fast speed, supporting instant printing.

Cassette Laser Printer KF-UE-006

PRODUCT FEATURES



Application Circumstance

Batch Printing and Single Quick Printing.



Strong Software Function

Simple to operate, drag-and-drop control interface, real-time preview of printed content, compatible with LIS, HIS, and PIS.



High-Stability System

Adopts LINUX system, safe and reliable, high-speed operation, support USB, network port, and WIFI interface connection.



High-Speed Printing

15/min



Contact Us

All the information in this material is for reference only, the company reserves all the rights, including updating and revision, sales representative shall be consulted on issues of parameter updating.

• Tel/WhatsApp/Wechat: +86 186-7233-4799

• Email: info@kfbio.cn

• Web: kfbio pathology.com

Technical Data

Cassette Laser Printer KF-UE-006



Laser printing
easy to maintain



Upload Slot Capacity
reaches 450 cassettes



Support offline independent
printing and network printing



Support linkage printing of
pathology information system

Cassette Type	Suitable for covered and uncovered cassettes, no need to specify laser type consumables
Upload Slot Capacity	6*75
Upload Method	Automatic rotary loading with 6 stations, 6 kind of colors can be printed at the same time
Download Slot Capacity	Continuously collect over 120 cassettes
Print Accuracy	2500dpi
Print Speed	15/min
Operation System	LINUX system
Operation Interface	≥10 inch, built-in colorful touch screen, easy to operate
System Interface	Support HIS, LIS, and PIS to realize data-sharing
Connection Type	Network port, USB, no need to connect computer, iPad or phone
Module Edit	Users can modify the print template content
Automatic Printing	Add data and each interface obtains data to realize automatic printing
Real-time Preview	Support real-time information display and cancel or pause printing
Records	Support data track to check previous contents, time, and operators
Dedusting Method	Equipped with professional dedusting device and replaceable filter cartridge to realize dust and odor elimination
Operating Voltage	100-240VAC, 50/60HZ
Printing Quality	With professional third-party quality and durability test report

